

Prometaphase Chromosome Preparation from Mouse Spleen (C57Bl/6)

Section of Cancer Genomics, Genetics Branch, NCI
National Institutes of Health

Reagents

Acetic acid, glacial

Antibiotic-Antimycotic 100x

10,000 U/ml Penicillin G sodium, 10,000 µg/ml streptomycin sulfate, 25 µg/ml amphotericin B

Gibco BRL, Cat. 15240-013

Colcemid, KaryoMAX Colcemid Solution, 10 µg/ml

Gibco BRL, Cat. 15210-016

Concanavalin A (5 µg/µl)

Sigma, Cat. C-5275

Fetal Bovine Serum (FBS) heat inactivated

Gibco BRL, 16140-022

L-Glutamine-200 mM, 100x

Gibco BRL, 25030-016

Homogenizer

Thomas Scientific, Cat. 3431D7

Lipopolysaccharides (LPS) 5mg

Sigma, Cat. L-2637

Methanol

Methotrexate, 500 mg

Sigma, Cat. M 8407

Potassium chloride (KCl)

RPMI Medium 1640

Gibco BRL, Cat. 21870-050

Preparation

Reagents

Amount

Concanavalin A

5 mg

Sterile water

1 ml

For a stock solution of 5µg/µl

RPMI Complete Medium

Components	Amount
RPMI Medium 1640	440 ml
Antibiotic-Antimycotic, 100X	5 ml
L-Glutamine-200 mM, 100X	5 ml
Fetal Bovine Serum (FBS)	50 ml

Fixative

Prepare fresh: methanol/acetic acid (3:1)

Hypotonic Solution

0.075M KCl in distilled water

Lipopolysaccharides (LPS)

Lipopolysaccharides (LPS)	25 mg
Sterile water	1 ml

Use 1:1000 dilution for a final concentration of 25 µg/ml of culture

MTX stock

Make an initial stock of 10^{-3} M in H₂O and then dilute to 10^{-5} M

BrdU stock

1 mg/ml in distilled water

Procedure

1. Prepare tissue culture flasks. To one T75 flask, add:

<u>Components</u>	<u>Amount</u>
Prepared media	20 ml
Concanavalin A (5µg/µl)	30 µl
Lipopolysaccharides (LPS)	25 µl

2. Isolate spleen from mouse. Transport in RPMI 1640.
3. Place three spleens into a homogenizer with 3 ml of plain RPMI media. Grind well.
4. Transfer 0.5 ml to one T75 flask.
5. Incubate at 37°C for 24 hr. After 24 h add 200 µl of MTX stock (10^{-5} M) to 20 ml of culture (MTX final concentration of 10^{-7} M); mix well and incubate an additional 17 h.

6. After 17 h centrifuge the content of the flasks, remove the supernatant, and wash the pellet twice with plain media.
7. After the second wash resuspend the pellet in 20 ml of RPMI 1640 10% BSA and transfer to a T75 flask.
8. Add 500 μ l of the BrdU stock (1mg/ml) to a final concentration of 25 μ g/ml (minimize light exposure).
9. Incubate for 5 h 30 min at 37°C.
10. For the last 10 min of the incubation add 20 μ l of Colcemid stock (10 μ g/ml) to a final concentration of 0.06 μ g/ml .
11. Centrifuge cultures for 10 min.
12. Transfer to 50 ml centrifuge tubes and centrifuge at 1000 rpm for 10 min.
13. Remove supernatant.
14. Gently add 10 ml 0.075M KCl (prewarmed to 37°C) to each tube and resuspend pellet.
15. Incubate tubes at 37°C for 15 min.
16. Following incubation, add a few drops of freshly prepared fixative.
17. Centrifuge at 1000 rpm for 10 min.
18. Remove supernatant.
19. Wash pellet with freshly prepared fixative, at least 3 times.
20. Store pellet under fixative at -20°C until ready to prepare slides.