

DNA Precipitation and Hybridization (CGH)

Section of Cancer Genomics, Genetics Branch, NCI
National Institutes of Health

Reagents

Dextran sulfate (50%)

Intergen, Cat. S4030

Ethanol, absolute**Formamide, deionized****Human Cot-1 TM DNA, 1 mg/ml**

GIBCO BRL, Cat. 15279-011, 500 µg

Salmon testes DNA, 9.7 mg/ml

SIGMA Molecular Biology, Cat. D-7657, 1 ml

SSC, 20X**Sodium acetate (Na-Acetate), 3M****Water, sterile**

Preparation

Master Mix

Dextran sulfate, 50%	40 ml
----------------------	-------

20X SSC, pH 7.0	10 ml
-----------------	-------

Sterile H ₂ O	50 ml
--------------------------	-------

Vortex solution and place tube on a shaking platform overnight to insure proper mixing.

Aliquot, and store at -20°C

70% Formamide/2X SSC

Deionized formamide	70 µl
---------------------	-------

20X SSC	3 µl
---------	------

Sterile water	27 µl
---------------	-------

Adjust to pH 7.5

Procedure

1. Add to an eppendorf tube:
 - 10-25 μ l nick-translated test or tumor probe DNA (200-500 ng DNA)
 - Equal amount of nick-translated control whole genomic DNA as probe DNA
(Note: can use 1-2 μ g DNA if tumor DNA is isolated from paraffin material)
 - 10-20 μ l human Cot-1 DNA (1 mg/ml)
 - 8-10 μ l salmon sperm DNA (10 mg/ml)
- Note:** Usually the test DNA is nick translated with Biotin-16-dUTP and the control DNA is nick translated with Digoxigenin-11-dUTP.
2. Add 1/10 volume Na-Acetate (3M).
3. Add 2.5-3.0 x total volume of absolute ethanol.
4. Vortex, store at -20°C overnight, or at -80°C for at least 15-30 min.
5. Centrifuge (14000 rpm) precipitated DNA at 4°C for 30 min.
6. Pour off supernatant and speed vac for 5-10 min to dry pellet.
7. Add 6 μ l deionized formamide (pH 7.5), incubate at 37°C for 30 min, shaking; vortex a few times during the 30 min incubation.
8. Add 6 μ l Master Mix, vortex, and centrifuge briefly.
9. Denature probe DNA at 76°C for 5 min and centrifuge briefly. Can keep at 37°C until ready to denature.
10. Preanneal at 37°C for 1-2 hr.
11. For slide denaturation apply 120 μ l 70% formamide/2X SSC to a 24 x 60 mm coverslip and touch slide to coverslip.
12. Incubate slides at 75°C for 1.5 min.
13. Immediately place in ice cold 70% ethanol for 3 min, followed by 90% ethanol and 100% ethanol for 3 min each; air dry.
14. Add probe DNA after preannealing to denaturated slides and cover with 18 mm² or 22 mm² coverslips; seal coverslips with rubber cement.
15. Hybridize at 37°C in a humidified chamber for at least 48 hr.

