

Immunocytochemistry Followed by FISH (Version 3)

Section of Cancer Genomics, Genetics Branch, NCI
National Institutes of Health

***We present multiple versions of the Immunocytochemistry / FISH protocols because the conditions under which each antibody and DNA probe work can be different and must be determined empirically.**

Reagents

Acetic acid, glacial

Bovine Serum Albumin (BSA)

Boehringer Mannheim Biochemicals (BMB), Cat. 100 350

DAPI

BMB, Cat. 236 276

Formamide

FLUKA BioChemica, Cat. 47670

Goat anti-rabbit-TRITC (secondary)

Sigma, Cat. T-5268

HCl, 1 M

Methanol

Para-Formaldehyde

Sigma, Cat. P6148

Phosphate Buffered Saline, pH 7.4

Invitrogen Corp., Cat. 10010-023

Primary antibody

Specific for desired protein, made in either a mouse or rabbit

Rabbit anti-mouse-TRITC (secondary)

Sigma, Cat. T2402

NaOH, 0.1 M

20X SSC

Tween 20

Sigma, Cat. P1379

Vysis CEP® Probe

Vysis

Preparation

Methanol

1. Room temperature
2. Pre-chill to -20°C

Permeabilization Buffer

Triton X-100	50 μ l	f.c. [0.53]
1X PBS	10 ml	

Blocking Solution (3% BSA/1X PBS)

BSA	0.3 g
1X PBS	10 ml

Store at 4°C

Antibody Solution (1% BSA/1X PBS)

Blocking solution	300 μ l
1X PBS	600 μ l

2% p-formaldehyde

p-formaldehyde	2 g	
1X PBS	100 ml	
0.1 N NaOH	500 μ l	f.c. [0.5 mM]
Adjust to pH 7.4 with HCl		
Store <1 month at 4°C		

50% FA/SSC

20X SSC	20 ml
dH ₂ O	80 ml
Formamide	100 ml
Total	200 ml

Adjust pH to 7-7.5 with 1 M HCl

Pre-warm to 45°C

DAPI (stock solution)

DAPI	2 mg	f.c. [0.2 mg/ml]
dH ₂ O	10 ml	
Aliquot and store at -80°C		

DAPI (staining solution)

DAPI stock solution	40 µl	f.c. [80 ng/ml]
2X SSC	100 ml	

Antifade (1,4-phenylene-diamine)

See Antifade preparation procedure in CGH Protocols

Procedure

1. Grow adherent cells in chamber slides.
2. Fix cells in methanol pre-chilled to -20°C for 10 min at RT.
3. Wash 3 x 5 min 1X PBS at RT.
4. Permeabilize cells with 0.5% Triton X-100/PBS 5 min at RT.
5. Wash 3 x 5 min 1X PBS at RT.
6. Remove chamber and block slides with 120 µl blocking solution in hybridization chamber 30 min at 37°C.
7. Incubate with 1° Ab (rabbit or mouse) in 120 µl antibody solution in hybridization chamber at 37°C for 45 min.
8. Wash 3 x 5 min with 1X PBS at RT.
9. Incubate with 2° Ab [goat anti-rabbit-TRITC (1:200) or Rabbit anti-mouse-TRITC (1:200), respectively, in 120µl antibody solution] in hybridization chamber at 37°C for 60 min.
10. Wash 3 x 5 min 1X PBS at RT.
11. Fix with methanol:acetic acid (3:1) at RT 10 min.
12. 2% p-formaldehyde at RT for 1 min.
13. 70%, 90%, 100% ethanol series (3 min each).

Note: Can counterstain with DAPI at this point and mount slides with antifade to have a look at them and determine if Ab detection worked. Wash coverslips 3 x 5 min in 2X SSC before continuing with procedure.

14. Combine 1 μ l Vysis CEP® probe, 1 μ l water, and 7 μ l Vysis Hyb Buffer.
15. Add probe cocktail to slide, coverslip, and seal with rubber cement.
16. Denature slide 75°C for 5 min on slide warmer.
17. Incubate in hybridization chamber at 37°C overnight.
18. Remove rubber cement.
19. Wash in FA/SSC pre-warmed to 45°C for 21 min, shaking.
20. Stain for 2 min with DAPI.
21. Wash in 2X SSC for 10 min, shaking.
22. Mount with antifade.