

## **Detection (SKY)**

Section of Cancer Genomics, Genetics Branch, NCI  
National Institutes of Health

### **Reagents**

**Avidin-Cy5 (1.8 mg/ml)**

Jackson Immuno Research Lab, Cat. 003-170-083

**Bovine Serum Albumin (BSA)**

Roche Diagnostics, Cat. 100-350

**DAPI**

Sigma, Cat. 18860

**Ethyl alcohol, anhydrous**

**Formamide**

Fluka BioChemika, Cat. 47671

**HCl, 1 N**

**Mouse anti-digoxigenin (0.1 mg/ml)**

Sigma, Cat. D 8156

**Sheep anti-mouse Cy5.5 (1.0mg/ml)**

Rockland Immunochemicals, Cat. 610-113-121

**SSC, 20X**

**Tween 20**

Sigma, Cat. P-1379

**dH<sub>2</sub>O**

**Antifade (1,4-phenylene-diamine)**

Sigma, Cat. P1519, 100 g

### **Preparation**

**50% FA/2X SSC**

20X SSC        20 ml

dH<sub>2</sub>O         80 ml

Formamide    100 ml

Adjust pH to 7.25 using 1 N HCl

**Pre-warm to 45° C**

**1X SSC**

20X SSC        25 ml

dH<sub>2</sub>O         475 ml

**Pre-warm to 45° C**

**4X SSC/0.1%Tween 20**

20X SSC        200 ml  
dH<sub>2</sub>O        799 ml  
Tween 20        1 ml

**Pre-warm to 45° C**

**Blocking Solution (3% BSA/4X SSC/0.1% Tween 20)**

BSA                    0.3 g  
4X SSC/0.1%Tween 20        10 ml  
Vortex until dissolved

**Pre-warm to 37° C**

**Antibody Solution (1% BSA/4X SSC/0.1% Tween 20)**

0.1 g BSA  
10 ml 4X SSC/0.1%Tween 20

**Pre-warm to 37° C**

(or use 4X SSC/ 0.1% Tween 20; see note 4)

**DAPI stock solution (f.c.= 0.2 mg/ml)**

DAPI                2 mg  
dH<sub>2</sub>O                10 ml  
Aliquot and store at -80°C

**DAPI staining solution (f.c.= 80 ng/ml)**

DAPI (stock solution)        40 µl  
2X SSC                    100 ml  
Store at 4°C in a light-tight coplin jar.

**Procedure**

1. Carefully remove rubber cement surrounding coverslips with forceps. Pre-soak slide in formamide/2X SSC if rubber cement is difficult to remove.
2. Wash slides in 50% formamide/2X SSC for 3 x 5 min each, shaking, preferably in 45°C water-bath.
3. Wash slides in 1X SSC for 3 x 5 min, shaking.
4. Dip slides in 4X SSC/0.1% Tween 20; do not let them dry.
5. Add 120 µl of Blocking Solution (3% BSA/4X SSC/0.1%Tween20) to a 24 mm x 60 mm coverslip and incubate in a moist hybridization chamber at 37°C for 30 min.
6. Wash slides in 4X SSC/0.1% Tween 20 to wash off blocking solution, 5 min, shaking.

7. Spin all fluorescent dyes for 1 min at 13,000 rpm.
8. Combine the two antibodies, mouse anti-Dig and Avidin-Cy5, into the same eppendorf tube, and apply 120  $\mu$ l of antibody solution to a 24 mm x 60 mm coverslip. Each antibody should be diluted 1:200 in 1% BSA/4X SSC/0.1% Tween 20 or 4X SSC/0.1% Tween 20 (see note 4). Invert the slide onto the solution. Incubate the slides in a moist hybridization chamber at 37°C for 45-60 min.
9. Wash slides in 4X SSC/0.1% Tween 20, 3 x 5 min, shaking.
10. Add 120  $\mu$ l of the antibody (sheep anti-mouse Cy5.5, diluted 1:200 in 1% BSA/4X SSC/0.1% Tween 20 or 4X SSC/0.1% Tween 20). Incubate slides in a moist hybridization chamber at 37°C for 45-60 min.
11. Wash slides in 4X SSC/0.1% Tween 20, 3 x 5 min, shaking.
12. Stain slides for 5 min in the DAPI staining solution in a light-protected coplin jar at RT.
13. Wash slides with 2X SSC 3-5 min.
14. Dehydrate slides in ethanol series of 70%, 90%, and 100% for 3 min each, air-dry slides.
15. Apply 35  $\mu$ l of antifade solution, cover each slide with a 24 mm x 60 mm coverslip, and store in a light-protected container at 4°C until slide is imaged.

## Notes

1. Exposure of slides to ambient light should be minimized during all procedures.
2. Carefully remove coverslips during all procedures to minimize scratches.
3. Do not let the slide dry out between washing steps.
4. BSA may contribute to non-specific background.