

Isolation of Interphase Nuclei from Paraffin Section for FISH, Hedley Technique

Reagents

EDTA, 0.5 M
Ethanol, absolute
Formaldehyde, 37%
Nytal HD 55 filter (70 µm)
1X Phosphate Buffered Saline (PBS)
Proteinase K
EM Science, Gibbstown, Cat. 24568-
NaCl solution, 0.9% and 5 M
Tris-HCl, pH 7.5
Xylene
Water, sterile

Preparation

| Proteinase K | |
|--------------------------|---------|
| Proteinase K | 0.005 g |
| sterile H ₂ O | 908 µl |
| Tris-HCl, pH 7.5 | 50 µl |
| 0.5 M EDTA | 20 µl |
| 5 M NaCl. | 2 µl |

| 1% Formaldehyde/PBS | |
|----------------------------|---------|
| Formaldehyde | 2.7 ml |
| 1X PBS | 97.3 ml |

Procedure

1. Remove excess wax by incubating slides at 65°C for 2 hr.
2. Transfer slides immediately to jar containing xylene. Incubate at room temperature 2 x 10 min.
3. Wash slides 5 min each in 100%, 86%, and 70% ethanol.

4. Wash slides 2 x 2 min in 0.9% NaCl.
5. Transfer sections to a petri dish and mince the tissue by cutting with scissors.
6. Transfer to eppendorf tube containing proteinase K solution, incubate for 30 min at 37°C.
7. Filter suspension using a Nyltal HD 55 filter (70 µm); wash the filter with 4 ml 1X PBS.
8. Centrifuge for 8 min at 1000 rpm.
9. Wash with 2 ml of 1X PBS, centrifuge for 8 min at 1000 rpm.
10. Resuspend pellet in about 20 µl of 1X PBS and drop onto clean glass slide or use a cytospin if possible.
11. Dry on a slide warmer at 37°C and then at room temperature overnight.
12. Rinse slides in PBS. Fix in 1% formaldehyde/PBS for 5 min, wash in PBS.
13. Drain excess fluid, dehydrate through ethanol series (70%, 90%, 100%), and air-dry.
14. Apply hybridization solution, add coverslip, and seal.
15. Denature the specimen at 85°C for 5 min and hybridize at 37°C for at least 16 hr. (Note: When probes are used that require a preannealing step, the denaturation of specimen and probe is performed separately.)