

Detection (FISH)

Introduction

Probes labeled with biotin must be detected with fluorescently labeled Avidin, and probes labeled with digoxigenin require detection with a fluorochrome conjugated antibody against this hapten. For example, to detect the biotin labeled probes we routinely use Avidin-FITC, Avidin-TRITC, or Avidin Cy-5. For the probes labeled with digoxigenin, we usually first incubate with mouse-anti-digoxigenin, followed by incubation with sheep anti-mouse Cy5.5, or other fluorochrome conjugated antibodies.

Reagents

Avidin-Cy5

Jackson Immuno Research Lab Cat 003-170-083

Avidin-TRITC

Sigma Cat. A 7169

Avidin-FITC

Vector Cat. A-2011

BSA (Bovine Serum Albumin)

DAPI

Ethanol, absolute

Formamide

Fluka BioChemika Cat 47671

HCl 1N

Mouse anti-digoxigenin

Sigma, Cat. D 8156

Sheep anti-mouse Cy5.5

Amersham, Cat. RPQ 0115

SSC 20X

Tween 20

Preparation of Reagents

50% FA/SSC

SSC 20X 30 ml

dH₂O 120 ml

Formamide 150 ml

Adjust pH to 7 with 1N HCl

Pre-warm to 45°C

1X SSC (for direct labeled probes, i.e., TRITC, FITC or other)

SSC 20X 25 ml

dH₂O 475 ml

Pre-warm to 45°C

0.1X SSC (for indirect labeled probes, i.e. Biotin, or Digoxigenin)

SSC 20X 2.5 ml

dH₂O 497.5 ml

Pre-warm to 60°C

4X SSC/0.1% Tween20

SSC 20X 200 ml

dH₂O 799 ml

Tween 20 1 ml

Pre-warm to 45°C

Blocking Solution (3% BSA/4X SSC/0.1% Tween20)

BSA 0.3 g

4X SSC/0.1% Tween 20 10 ml

Pre-warm to 37°C

Antibody Solution (1% BSA/4X SSC/0.1% Tween 20)

BSA 0.1 g

4X SSC/0.1% Tween 20 10 ml

Pre-warm to 37°C

DAPI stock solution (f.c.= 0.2mg/ml)

DAPI 2 mg

ddH₂O 10 ml

Aliquot and store at -80°C

DAPI staining solution (f.c.= 80ng/ml)

DAPI (stock solution) 40 µl

SSC 2X 100 ml

Store at 4°C in a light-tight coplin jar

Procedure

1. Carefully remove the rubber cement surrounding the coverslips from hybridized slides.
2. Wash the slides in 50% formamide/2xSSC (pH 7-7.5) for 3 x 5 min, shaking.
3. Wash slides in 0.1X SSC (for indirectly labeled probes) or 1X SSC (for directly labeled probes) for 3 x 5 min, shaking.

4. Dip slides in 4X SSC/0.1% Tween 20.
5. Add 120 μ l of Blocking Solution (3% BSA/4X SSC/0.1% Tween 20) to the slides and over them with a 24 mm x 60 mm coverslip in a moist hybridization chamber at 37°C for 30 min.
6. Dip slides in 4X SSC/0.15% Tween 20 to wash off the blocking solution.
Proceed directly to step 9 if using a directly-labeled probe.
7. For indirectly-labeled probes (Biotin or Digoxigenin), add 120 μ l of fluorescent antibody (antibody should be diluted 1:200 in 1% BSA/4X SSC/0.1% Tween 20) to the slides, cover with a 24 mm x 60 mm coverslip, and incubate in moist light-tight hybridization chamber at 37°C for 45 min.
8. Wash slides in 4X SSC/0.1% Tween 20, for 3 x 5 min, shaking.
9. Stain slides for 5 min in DAPI staining solution in a light-protected coplin jar.
10. Wash the slides for 5 min in 2X SSC, shaking.
11. Dehydrate the slides by dipping through 2X SSC followed by washes with an increasing ethanol series of 70%, 90%, and 100%; air-dry.
12. Apply 35 μ l of antifade solution, cover with 24 mm x 60 mm coverslips, store in light-protected container at 4°C until slide is imaged.

Notes

1. Exposure of slides to ambient light should be minimized during all procedures.
2. Use care in removing coverslips during all procedures to minimize scratches.
3. Spin all fluorescent dyes prior to use for 3 min at 13,000 rpm and carefully pipette the antibody without disturbing the pellet..
4. Do not let the slide dry out between washing steps.