

# Nick Translation (FISH)

## Reagents

**Biotin-16-dUTP (100mM) (see Notes)**

Boehringer Mannheim, Cat. 109 3070

**Bovine serum albumin (BSA)**

**dATP, dTTP, dGTP, dCTP**

Boehringer Mannheim, Cat.: 105 1440, 105 1458, 105 1466, 105 1482

**DNase I from bovine pancreas**

Boehringer Mannheim, Cat. 104 159, 100 mg

**99.5% Glycerol**

Gibco BRL, Cat. 15514-011

**Magnesium chloride (MgCl<sub>2</sub>) (0.5 M)**

**β-Mercaptoethanol Solution (99%)**

Sigma, Cat. M6250

**DNA Polymerase (Kornberg)**

Boehringer Mannheim, Cat. 104 485, 500 U

**NaCl (1 M)**

**Tris-HCl (1 M), pH 8.0**

Quality Biologicals, Cat. 351-007-100

**H<sub>2</sub>O**

## Preparation

**Bovine DNase stock solution, 1mg/ml**

Dnase, 10 mg

NaCl (1M), 1.5 ml to get a final concentration of 0.15 M

Glycerol (100%), 5 ml to get a final concentration of 50%

H<sub>2</sub>O to total volume of 10 ml

Mix well and store at –20°C.

**0.1M β-Mercaptoethanol**

– Mercaptoethanol solution (99%/14.4 M), 34.7 µl

H<sub>2</sub>O to a total volume of 5 ml

Mix well aliquot and store at –20°C.

**dNTP**

dATP (100 mM), 5  $\mu$ l for a final concentration of 0.5 mM

dCTP (100 mM), 5  $\mu$ l for a final concentration of 0.5 mM

dGTP (100 mM), 5  $\mu$ l for a final concentration of 0.5 mM

dTTP (100 mM), 1  $\mu$ l for a final concentration of 0.05 mM

H<sub>2</sub>O, 984  $\mu$ l

Mix well, aliquot, and store at  $-20^{\circ}\text{C}$

**10X NT-Buffer**

Tris-HCL (1 M, pH 8.0 ), 500  $\mu$ l for a final concentration [0.5 M]

MgCl<sub>2</sub>(0.5 M), 100  $\mu$ l for a final concentration [50 mM]

BSA (10 mg/ml), 50  $\mu$ l for a final concentration [0.5 mg/ml]

H<sub>2</sub>O, 350  $\mu$ l

Mix well, aliquot, and store at  $-20^{\circ}\text{C}$ .

**Procedure****1. Combine components for 100  $\mu$ l Nick Translation Reaction**

DNA	X $\mu$ l (2 $\mu$ g)
10X NT Buffer	10 $\mu$ l
0.1M $\beta$ -mercaptoethanol	10 $\mu$ l
dNTP	10 $\mu$ l
BIO-16-dUTP	4 $\mu$ l
Sterile H <sub>2</sub> O	60 $\mu$ l – X $\mu$ l
Diluted Dnase (note)	4 $\mu$ l
Polymerase (Kornberg)	2 $\mu$ l

Total volume of 100  $\mu$ l.

Mix and briefly centrifuge tube.

2. Place tube into a  $15^{\circ}\text{C}$  waterbath for 1 hr for the nick translation reaction.
3. Stop the reaction by placing tube on ice.
4. Run gel with 5  $\mu$ l of sample and 0.5 ml of loading buffer.
5. Length of DNA should be about 300 bp-600 bp. See image.

## Notes

1. Although Biotin is the only analog listed as the label for the DNA, other analogs can be used such as Digoxigenin-11-dUTP (Boehringer Mannheim, Cat. 155 8706)
2. DNase should be diluted 1:1000 for the initial use and adjusted up or down according to the reaction velocity. See the image. If the DNase reaction velocity is too high, it will result in fragments that are too small. See lanes 7,8,9. One should increase the dilution factor (for example 0.8:1000). If the DNase reaction velocity is too slow it will result in fragments that are too large. See lane 2. One should decrease the dilution (for example 1.2:1000).

